

## Rehabilitation Protocol: Meniscus Allograft Transplantation

Name:	Date:
Diagnosis:	Date of Surgery:
Phase I (Weeks 0-8)	
• Weightbearing:	
• Weeks 0-4: Toe touch v	weighthearing
o Weeks 4-6: Advance to	WBAT with crutches (d/c crutches at 4 weeks post-op if gait normalized)
Hinged Knee Brace: worn for 6	6 weeks post-on
<ul> <li>Locked in full extension</li> </ul>	for ambulation and sleeping – remove for hygiene (Week 1)
<ul> <li>Locked in full extension</li> </ul>	for ambulation- remove for hygiene and sleeping (Weeks 2-4)
<ul> <li>Set to range from 0-90°;</li> </ul>	for ambulation- remove for hygiene and sleeping (Weeks 4-6)
<ul> <li>Discontinue brace at 6 w</li> </ul>	veeks post-op.
<ul> <li>Range of Motion – PROM → AA</li> </ul>	AROM → AROM as tolerated
<ul> <li>Weeks 0-4: Non-weight</li> </ul>	tbearing 0-90°
<ul> <li>Weeks 5-8: Full non-we</li> </ul>	eightbearing ROM as tolerated – progress to flexion angles greater than 90°
<ul> <li>Therapeutic Exercises</li> </ul>	
<ul> <li>Quadriceps sets, heel slid</li> </ul>	des, straight leg raises, patellar mobilizations, co-contractions (Weeks 0-2)
<ul> <li>Add neel raises and term</li> </ul>	ninal knee extensions (Weeks 2-8)
<ul> <li>Activities in brace for fir</li> </ul>	st 6 weeks – then without brace
<ul> <li>No weightbearing with</li> </ul>	n flexion > 90° during weeks 0-4
<ul> <li>Avoid tibial rotation fo</li> </ul>	or first 8 weeks to protect the meniscal allograft
Phase II (Weeks 8-12)	
Weightbearing: As tolerated	
• Range of Motion – Full active R	ROM
Therapeutic Exercises	
	extension exercises, begin hamstring strengthening
o Lunges – 0-90°, Leg pres	s = 0.90° (flexion only)
o Proprioception exercises	s
o Begin use of the stationa	
Phase III (Months 3-6)	and the state of t
Weightbearing: Full weightbea     Range of Motion - Full / Painless	iring with normal gait pattern
ango or motion Tuny Lannes.	S ROM
Therapeutic Exercises	
o Continue with quad and	namstring strengthening
o Focus on single-leg stren	igth
o Begin jogging/running	
o Plyometrics and sport-sp	
Maintenance magazine for attrict	y as tolerated (6 months post-op)
<ul> <li>Maintenance program for streng</li> </ul>	gth and endurance
Comments:	
Frequency: times per week	Duration: weeks
o	
Signature:	Date: